

#655N/#657N Care Cliners



ISO 9001-2000 CERTIFIED

All of Our Best Features - Plus the Added Versatility of Drop Arms

Leave it to Winco to combine the best clinical chair features in the industry with the convenience of dual drop-arms!

Our Drop-Arm Care Cliners make maneuvering your patients easier and less stressful - for both you and them. Our removable side-panels are standard and make clean-ups quick and easy.



#657N
XL Droparm
Care Cliner

#655N
Dual Drop-Arm
Care Cliner

For Use in Dialysis,
Oncology, Acute Care,
& Other Clinical Environments



Proudly Made in the USA!

1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474

www.wincomfg.com

#655N/#657N Care Ciner



ISO 9001-2000 CERTIFIED

Standard Trendelenberg Release on Back with Easy-Access Handle

Built-in, Ergonomic Headrest

Webbed Elastic and an Open Seat Base Provide Superior Comfort and Support

Adjustable Elastic Bands on Back Frame Provide Exceptional Comfort

Non-Porous, Fold-Down Side Trays Contain Small Spills and are Standard on Both Sides

One-Touch Release for Drop-Arms

Removable Side-Panels for Easy Cleaning and Maintenance

High-Impact, Nylon, 5" (12.7cm) Casters - Directional-Lock in Left Rear. -All others are Total-Lock.

VINYL SELECTION

(Printed vinyl colors may vary slightly from actual)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
02	03	04	06	07	14	16	17	OPTIONAL 32	OPTIONAL 33
OPTIONAL 34	OPTIONAL 37	OPTIONAL 41	OPTIONAL 42	OPTIONAL 43	OPTIONAL 44	OPTIONAL 60	OPTIONAL 61	OPTIONAL 62	OPTIONAL 63
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress


#655N SPECIFICATIONS

#657N SPECIFICATIONS

OPTIONS

SEAT HEIGHT: 21" (53.3cm)
 SEAT WIDTH: 21" (53.3cm)
 SEAT DEPTH: 20" (50.8cm)
 OVERALL HEIGHT: 49" (124.5cm)
 OVERALL WIDTH: 33" (83.8cm)*
 BACK HEIGHT: 31.5" (80cm)
 WEIGHT: 105 lbs. (47.6 kg)
 WEIGHT CAPACITY: 275 lbs. (124.7 kg)
 FOOTPRINT: 33" x 40" (83.8cm x 101.6cm)*
 FOOTPRINT (RECLINED): 33" x 76.5" (83.8cm x 194.3cm)*

SEAT HEIGHT: 21" (53.3cm)
 SEAT WIDTH: 25" (63.5cm)
 SEAT DEPTH: 20" (50.8cm)
 OVERALL HEIGHT: 49" (124.5cm)
 OVERALL WIDTH: 36" (91.4 cm)
 BACK HEIGHT: 31.5" (80cm)
 WEIGHT: 111 lbs. (50.4 kg)
 WEIGHT CAPACITY: 450 lbs. (204 kg)
 FOOTPRINT: 36" x 43" (91.4cm x 109.2cm)
 FOOTPRINT (RECLINED): 36" x 76.5" (91.4cm x 194.3cm)



TB-133 Where Required
#655NT #657NT



Steel Casters
#655 #657



Heat & Massage
#655NH #657NH
 Heat ONLY
#655NW #657NW



IV Pole
#IH #IH



Headrest Cover
#HC #HX

* ADD 9" (22.9cm) TO WIDTH PER SIDE TABLE WHEN RAISED.