

# #656 Premier Care Cliner



ISO 9001-2000 CERTIFIED

*Style, Comfort, Value, and Great Features Make  
Our Premier Care Cliner Unlike any Other*

Our Premier Care Cliner has a uniquely smooth, contemporary design to help your patients feel at ease, while maintaining the quality, durability, and great features you've come to expect from a Winco. It's no wonder Winco Care Cliners are among the best selling clinical chairs on the market today!



#656  
Premier  
Care Cliner

For Use in Dialysis,  
Oncology, Acute Care,  
& Other Clinical Environments



Proudly Made in the USA!

1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474


[www.wincomfg.com](http://www.wincomfg.com)

# #656 Premier Care Cliner



ISO 9001-2000 CERTIFIED

Adjustable Elastic Bands on Back Frame Provide Exceptional Comfort



Elastic Open Seat Base Provides Superior Comfort and Support

Non-Porous, Fold-Down Side Tray Contains Small Spills  
Standard on Right Side  
Optional on Left Side



True Trendelenberg Positioning Comes Standard

Low-Profile, Total-Lock, Steel Casters Provide Smooth, Easy Maneuvering

## VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
<b>02</b>	<b>03</b>	<b>04</b>	<b>06</b>	<b>07</b>	<b>14</b>	<b>16</b>	<b>17</b>	<b>32</b>	<b>33</b>
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress
<b>34</b>	<b>37</b>	<b>41</b>	<b>42</b>	<b>43</b>	<b>44</b>	<b>60</b>	<b>61</b>	<b>62</b>	<b>63</b>

## #656N SPECIFICATIONS

SEAT HEIGHT:	20" (50.8cm)
SEAT WIDTH:	21" (53.3cm)
SEAT DEPTH:	19.5" (49.5cm)
OVERALL HEIGHT:	46.5" (118.1cm)
OVERALL WIDTH:	30.5" (77.5cm)*
BACK HEIGHT:	30" (76.2cm)
WEIGHT:	93 lbs. (42.2 kg)
WEIGHT CAPACITY:	275 lbs. (124.7 kg)
FOOTPRINT:	30.5" x 43" (77.5cm x 109.2cm)*
FOOTPRINT (RECLINED):	30.5" x 76.5" (77.5cm x 194.3cm)*

\* ADD 10" (25.4cm) TO WIDTH PER SIDE TABLE WHEN RAISED.

## OPTIONS



TB-133  
Where Required  
**#656T**



2nd Side Table On Left  
**#656D**



IV Pole  
**#IM**



Headrest Cover  
**#HC**