


#570/#571 Blood-Drawing Chairs



ISO 9001-2000 CERTIFIED

Value and Quality at an Affordable Price

The high quality features of our Blood Drawing Chairs, such as durable seat and back, laminated cabinet and drawer with glides, and sturdy leveling feet, will provide you with years of trouble-free service.



#570
Blood Drawing Chair
w/ Cabinet



#571
Without Cabinet



Proudly Made in the USA!


1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474

www.wincomfg.com

#570/#571 Blood-Drawing Chairs



ISO 9001-2000 CERTIFIED



Optional IV Pole
(#570-IV or #571-IV)

Durable Plastic Resin
Seat Stands Up to
Years of Use

Optional Pivot Arm Available on both
Models - Order #570P or #571P



Fold-Up Armrest Provides
Easy Access and Added
Safety for the User

Fully Laminated Cabinet
and Storage Drawer

Adjustable Height
Armrest

Powder-Coated
Steel Frame

High Impact Leveling
Feet Ensure a Sturdy
Chair on Almost Any
Surface

New Higher
20.5" (52cm)
Seat Height

Reinforced Cross Brace Keeps
your Chair Strong and Stable

#570 SPECIFICATIONS

#571 SPECIFICATIONS

SEAT HEIGHT:	20.5" (52CM)
SEAT WIDTH:	18" (45.7CM)
SEAT DEPTH:	15" (38CM)
OVERALL HEIGHT:	36" (91.4CM)
OVERALL WIDTH:	35" (88.9CM)
BACK HEIGHT:	16" (40.6CM)
WEIGHT:	44LBS. (20KG)
WEIGHT CAPACITY:	300LBS. (136KG)
FOOTPRINT:	35" x 20.5" (88.9cm x 52cm)
ARM HEIGHT RANGE:	28" x 38" (71.1cm x 96.5cm))

SEAT HEIGHT:	20.5" (52CM)
SEAT WIDTH:	18" (45.7CM)
SEAT DEPTH:	15" (38CM)
OVERALL HEIGHT:	36" (91.4CM)
OVERALL WIDTH:	29" (73.7CM)
BACK HEIGHT:	16" (40.6CM)
WEIGHT:	30LBS. (13.6KG)
WEIGHT CAPACITY:	300LBS. (136KG)
FOOTPRINT:	29" x 20.5" (73.7CM x 52CM)
ARM HEIGHT RANGE:	28" x 38" (71.1CM x 96.5CM))