


#572/#573/#575
Padded Blood-Drawing Chairs




ISO 9001-2000 CERTIFIED

*All the Features you need,
All the Quality you Demand*

Features like our contoured, adjustable, fold-up armrest, convenient side cabinet, and comfortable seat and back make our blood-drawing chairs your first choice. Our outstanding quality will provide you with years of trouble-free service.



#572
Padded
Blood-Drawing Chair
w/ Cabinet



#573
Without
Cabinet



#575
XL Chair w/o Cabinet



Proudly Made in the USA!

1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474

www.wincomfg.com

#572/#573/#575 Padded Blood-Drawing Chairs



ISO 9001-2000 CERTIFIED

Optional Pivot Arm Available



Fully Laminated Storage
Drawer with Padded Top

Comfortable
Padded Seat and
Backrest

New Higher
22.5" (57.2cm)
Seat Height

Durable, Powder-Coated
Steel Frame

Reinforced Cross Brace Keeps
your Chair Strong and Stable

Fold-Up Padded Armrest
Provides Easy Access and
Added Safety for the User

Adjustable Height
Armrest

High Impact Leveling
Feet Allow the Chair to
Remain Level on Most Surfaces

VINYL SELECTION

(printed vinyl colors may vary slightly from actual colors)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
02	03	04	06	07	14	16	17	OPTIONAL 32	OPTIONAL 33
OPTIONAL 34	OPTIONAL 37	OPTIONAL 41	OPTIONAL 42	OPTIONAL 43	OPTIONAL 44	OPTIONAL 60	OPTIONAL 61	OPTIONAL 62	OPTIONAL 63
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

572/#573/#575 SPECIFICATIONS

OPTIONS

#572

#573

#575

SEAT HEIGHT:	22.5" (57.2CM)	22.5" (57.2CM)	22.5" (57.2CM)
SEAT WIDTH:	18.75" (47.6CM)	18.75" (47.6CM)	24.25" (61.6CM)
OVERALL HEIGHT:	44.5" (113CM)	44.5" (113CM)	36.5" (92.7CM)
ARMREST HEIGHT:	28.5"-38.5" (72.4CM-97.8CM)	28.5"-38.5" (72.4CM-97.8CM)	28.5"-38.5" (72.4CM-97.8CM)
OVERALL WIDTH:	35" (88.9CM)	29" (73.7CM)	34" (86.4CM)
WEIGHT:	65 LBS. (29.5 KG)	48 LBS. (21.8KG)	53 LBS. (24KG)
WEIGHT CAPACITY:	300 LBS. (136 KG)	300 LBS. (136KG)	450 LBS. (204KG)
FOOTPRINT:	35" X 27" (88.9CM X 68.6CM)	29" X 27" (73.7CM X 68.6CM)	34" X 25" (86.4CM X 63.5CM)



IV Option
#572IV
#573IV
#575IV



Pivot Arm
#572P
#573P
#575P