


Standard AmeriGlide Bath Lifts are smaller units than Premium Bath Lifts but still have the "heavier than it looks" feel. This feel is the most significant sign of the bath lift's high quality construction. It's a one piece bathtub lift, and all you have to do to start using it is place it in your tub and connect the controller. The controller is identical to the one used on the AmeriGlide Premium bathtub lift and has large buttons that make it a joy to use.

With bath lifts, you can easily care for a family member in your home. Or, you can continue your privacy and independence without needing to ask for assistance. Enjoy the added safety with a bathtub lift from AmeriGlide.

AmeriGlide Bath Lift Specifications

Lifting Capacity:	300 lbs
Power Source:	Battery
Seat Width:	14"
Seat Width - Flaps Down:	28.5"
Seat Depth:	16.1"
Seat Height, Lowered:	2.8"
Seat Height, Raised:	18.3"
Backrest Width:	14"
Backrest Height:	25.8"
Weight, Unit:	24 lbs
Rotating Seat:	Optional
Installation Method:	Suction Cups
Warranty:	One Year Complete
Reclines:	No
Padding:	Optional back and seat covers

AmeriGlide Bath Lift Accessories



BL105 Swivel Slide

Eases entry and exit from your bath. Our most popular accessory.




BL110 Chest Harness with Pommel

Ensures occupant safety by preventing any slippage for patients with compromised upper body strength.



BL115 Lap Harness with Pommel

Ensures occupant safety by preventing downward slippage.



BL120 Headrest with Wings

Provides extra support and security for the head and upper body.



AmeriGlide Bath Lift Key Features

- Comfortable, washable padding
- Battery operated with low-battery safety switch
- Suction cup feet for stability
- Waterproof, floating hand controller with large buttons
- Simple to disassemble for storage and transportation
- 300lb. weight capacity
- Easy to use battery charger included
- Many accessories available
- Great price