

BACKPACKER ™ INTERIOR LIFT

Your on-the-go lifestyle demands that you bring your scooter or power chair with you. Now, Pride Mobility Products® Corp., the number one power chair and scooter manufacturer, offers the lift you need to take your mobility product on the road. Backed by the service you have come to expect and trust from Pride®, with the full support of our expert technical service staff, an AmeriGlide™ Lift is the natural complement to your Pride Scooter or Power Chair.

Take your mobility on the road with the AmeriGlide Backpacker[™]interior lift from Pride. It's the ultimate solution for lifting your power chair or scooter into your full size van or minivan for transport. Like all AmeriGlide Lifts from Pride, the Backpacker features all-steel construction and a protective finish, along with an industrial rated 12 volt motor and optional battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

A Division of **Pride** Mobility Products Corp.

BACKPACKER ™ INTERIOR LIFT







Backpacker fully loaded.

Backpacker in lowest position.

Backpacker raising.

Base weight lift	106 lbs. ¹
Scooter / power chair overall weight	325 lbs.
Maximum product width	26"
Maximum product wheelbase	35"
Maximum travel distance from the ground	26"
Motor	Sealed 12 volt DC
Warranty	3-year transferable

Can easily lift the following Jazzy® Power Chairs ²

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1107
- CIFICATION Jazzy® 1113 ATS
 - Jazzy® 1121
 - Jazzy® 1122
 - Jazzy® Select
 - Jazzy® Select 14
 - Jazzy® Select 14XL

Can easily lift the following Quantum® Power Chairs ²

Optional battery pack12 volt, on-board

• Quantum® 600

Optional charger110 volt

- Quantum® 610
- Quantum® 1107
- Quantum® 1121
- Ouantum® 1402

Can easily lift the following Pride® Scooters ²

- Celebrity® X 3 & 4
- Go-Go® Ultra 3 & 4
- Go-Go Elite Traveller™ 3 & 4
- Legend® 3 & 4
- Maxima 3 & 4
- Rally® 3 & 4
- Revo® 3 & 4
- Sonic®
- Sundancer®
- Victory® 3 & 4

Can easily lift the following Jet™ Power Chairs ²

- Jet™ 3 Ultra
- Jet[™] 7
- Jet[™] 10 Ultra
- · Go-Chair®
- Weight does not include any accessories
- Vehicle dependent

Note: The Backpacker TM can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. Can fit in most mini-vans with little or no vehicle modification, full-size vans require modification.

