

VICTORY 10

THE ULTIMATE 10" SCOOTER
WITH LUXURY FEATURES
AT AN ECONOMY PRICE

VICTORY® 10

The sleek, sporty **Victory® 10** scooter delivers high performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features like exclusive low-profile tires, LED battery gauge, and wraparound delta tiller, the Victory® 10 is a greater value than ever before.



Delta tiller with wraparound handles



Lightweight seat with foam inserts and sliders



Pride's exclusive stylish, light weight, non-marking low profile wheels



New LED battery meter for increased visibility



Long-lasting, high-intensity LED headlight offers optimal pathway illumination

#1 in Performance


#1 in Reliability

#1 in Value



Victory® 10	3-wheel	4-wheel
Model number:	SC610	SC710
Weight capacity:	350 lbs.	350 lbs.
Maximum speed ¹ :	Up to 5.25 mph	Up to 5.25 mph
Ground clearance ² :	3"	3"
Turning radius:	43.25"	54"
Overall length ² :	46"	47"
Overall width ² :	22.25"	22.25"
Range per charge ¹ :	Up to 25 miles	Up to 25 miles
Tire Size:	Front: 3.5x10 solid Rear: 3.5x10 solid	Front: 3.5x10 solid Rear: 3.5x10 solid
Total Weight Without Batteries:	130.5 lbs.	134.5 lbs.
Heaviest piece when disassembled:	Front: 54.5 lbs.	Front: 61 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width 17" depth	
Seat-to-Deck Range:	16.75" - 18.25"	
Drive System:	Rear-wheel drive, sealed transaxle, 24-volt DC motor	
Dual Braking System:	Electronic, regenerative, and electromechanical	
Battery Requirements ³ :	Two 12-volt, deep cycle Size: U-1 or optional 40 AH Weight: U-1:24.5 lbs. each; 40 AH: 32.5 lbs. each	
Battery Charger:	3-amp, o -board with U-1 batteries 5-amp, o -board with optional 40 AH batteries	

Colors:



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Features

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Pride's exclusive stylish, light weight, non-marking, low profile wheels
- New LED battery meter for increased accuracy and visibility in any setting
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight, seat with viscoelastic foam inserts for increased comfort and sliders
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Pride Silver Star Lifts & Ramps

- Backpacker® Plus
- Scooter double entry
- Full Platform



Backpacker® Plus

These lifts are compatible with the Victory® 10.

Pride[®]
Mobility Products Corp.