

# The Key Features




**Easy to use Joystick**  
Makes the use of the stairlift easier and more comfortable even for those with limited mobility.

Exclusive to Platinum



**Isolating Key**  
located under the armrest.



**Footrest lever**  
to lift and lower the footrest without having to bend down.



**Ergonomically designed seat**  
Designed for optimised comfort for the position of the legs whilst seated.

Exclusive to Platinum



**Seat covers removable and washable**  
to keep your stairlift clean and tidy.

Exclusive to Platinum

**InVISION display**  
set into the armrest to have all information at your fingertips.

Exclusive to Platinum

**Safety harness**  
for the maximum peace of mind.

**Swivel seat**  
reaching the top of the stairs you can alight straight onto the landing.

**Emergency hand winder**  
in an emergency the seat can be moved using a manual handle.

**ON/OFF switch**  
to switch the stairlift on and off as required.

**12 Security sensors**  
A number of safety sensor will stop the chair safely if an obstruction is detected.

