

UpEasy Power Seat

Product SKU Number: UPEP100

Description: This portable electric lifting seat helps you sit down and rise easily without exerting undue stress on joints and muscles. Simply place the lifting seat on almost any armchair or sofa, sit and press the lever to slowly lower yourself. To get up, press the lever again and the lifting seat curves to raise you on a level plane slowly and gently without dumping you forward. Supports 100% of your weight to help you stand with little effort and the lift may be stopped at any height. The integrated handle ensures portability and the non-skid bottom keeps the lifting seat stable. The stain-resistant polyurethane foam cushion has a washable cover. Supports up to 300 lbs. Plugs into AC. 19" L x 17" W x 4" H. (12 lbs.)

Product Features

- 1. CSA & C/US Approved: Power Seat CSA & C/US approved. 10ft (3m) power cord plugs into a standard 120 volt outlet.
- 2. Portable: Power Seat is easy to move from room to room, or take with you when you are visiting. The base is designed with a built-in carry handle.
- 3. Durable: The injection-molded base is made of durable polycarbonate, and covered with a non-slip material that prevents shifting on fabric seating surfaces.
- 4. Washable: The marine blue fabric cover is made of 100% polyester and is machine washable in cold water.
- 5. Easy Care: The foam cushion has a water resistant coating that can be wiped down with mild soap and water.
- 6. Safe: The flexible seat curves to a flat surface that lifts you up on a level plane. This patented LeveLiftT design provides for a safe lift up, with no forward dumping motion. Patent # 5,316,370
- 7. Cushioning: The high density, polyurethane foam cushion is molded to fit body contours.
- 8. Easy to Use: Simply click the power lever up to rise, push down to lower. The middle position lets you stop at any time or height you desire.
- 9. Power Levers: Power Seat comes with two detachable power levers that can be used on either the left or right side on the cushion, according to your preference or need.



Pictures of Product:





Product Classification: Mobility Aids, Lifting Seats