

Specifications

o p c c i i c a c i c	119
Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.5"
Turning radius	24.75"
Overall length	42" without front riggings
Overall width	23"
Drive wheels	9" solid
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 8.95 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight Base weight Standard Seat weight	18 lbs. each Front: 47.6 lbs. Rear: 21 lbs. 43 lbs.
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of \pm or \pm %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Features

- · Easy disassembly for convenient transport
- · In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- · Simple main frame design for easy serviceability

Options

- High-back seat with solid seat pan
- Depth-adjustable, limited recline high-back
- · Cup holder
- · Weather cover
- · Rear basket
- · Oxygen tank holder
- Cane/crutch holder
- · Walker holder

Recommended Lifts

- Backpacker[®] Plus
- Outlander
- Backpacker[®] MV

Color Choices





(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com