




**Pride® Lift Chairs**  
*Experience the Luxury*

Protect your investment with the Pride® Lift Chair Care Kit



The Classic Collection from Pride® Lift Chairs delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection Lift Chairs feature patented, quiet and smooth operation and stylish designs to add to any home décor.

# *Classic* COLLECTION

[www.pridemobility.com](http://www.pridemobility.com)

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

**Pride®**  
Mobility Products Corp.

**CLASSIC COLLECTION MODELS**



CL-10	CL-15 <sup>2</sup>	CL-20	CL-30 <sup>2</sup>
2-Position Medium Button-back Design 5'2" - 5'10"	3-Position Medium Button-back Design 5'2" - 5'10"	2-Position Medium Sewn Pillow-back Design 5'4" - 6'0"	3-Position Medium Chaise Lounger Split-back Design 5'4" - 6'0"

**SPECIFICATIONS**

	CL-10	CL-15	CL-20	CL-30
<b>Seat-To-Floor:</b>	20"	20"	20"	20"
<b>Seat Depth:</b>	20"	20"	20"	20"
<b>Seat Width:</b>	20"	20"	20"	20"
<b>Top of Back to Seat:</b>	25"	27"	24.5"	26"
<b>Distance From Wall:</b>	12"	18"	11"	18"
<b>Overall Chair Width:</b>	33.5"	33.5"	32.5"	34"
<b>Chair Weight:</b>	110 lbs.	119 lbs.	114 lbs.	118 lbs.

**BENEFITS**

- Engineered furniture grade laminate/ hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head & arm covers
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets
- Removable back - no tools required
- 375 lbs. weight capacity on CL-20 & CL-30; 325 lbs. weight capacity on CL-10 & CL-15
- Integrated battery backup in case of power failure<sup>1</sup>
- Transformer features self-diagnostic electronics with 12' cord
- Lifetime warranty on steel lift frame
- 2-year limited warranty

**FABRIC SELECTION**

**STANDARD FABRICS**

MARINE WINE MOSS WHEAT TAUPE

---

**OPTIONAL FABRICS**

NEW RASPBERRY HERB NEW ALMOND TRUFFLE RASPBERRY BRASS ROSE TWIST NEW OATMEAL SLATE

**PADDED-SUEDE** BLUE EXPRESSO **PLUSH VELVET COLLECTION** BERRY CHOCOLATE **WINDSOR MICRO-SUEDE COLLECTION** CAMEL MAPLE

**LEXIS VINYL** BLACK MUSHROOM BURGUNDY NEW CHESTNUT ROYAL BLUE **NEW OUTLAST\* TECHNOLOGY** CHAMPAGNE NUTMEG

**CRYPTON® SUPER FABRICS SMART SUEDE\*\*** OYSTER PICKLE CLOVER MALT GINGER GUMBALL SOY TERRA COTTA CAPPUCCINO

**SMART FABRICS\*\*** VANILLA BEAN RED PEPPER CELERY CREAM

**CRYPTON Super Fabrics**

**OPTIONS\***

- Heat & Massage (CL-20 & CL-30 only)
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers
- Leather/vinyl wipes
- Footrest extension available on CL-15 & CL-30
- Crypton Malt & Cappuccino also available as CAL133 on CL-30

Additional lead time may apply for some options.

\* Additional charges apply for all options.

1. Batteries not included.
2. Available with optional footrest extension.

Look for the Quick Ship Indicator for the shortest lead times available.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

\*\*5-year warranty.