

Pride® Lifts

Outlander™

The **Outlander™** Exterior Lift was specifically designed with travel mobility scooters in mind for an easy-to-use power operated lift to take your travel mobility scooter with you. The open platform reduces the overall weight of the lift system to increase compatibility with Class II hitches. The lift features an integrated automatic scooter lockdown arm designed to provide a secure self-locking mechanism for compatible mobility scooters. Simple, hitch-based installation requires little to no modification of your vehicle and allows for easy removal.

FEATURES

- 150 lb (68 kg) lifting capacity
- Compatible with Class II, III, and IV hitches
- Compatible with Travel Mobility Scooters
- LED puddle light for ease of use at night
- Phosphate undercoating for weather resistance
- Manual override with tool kit
- Standard hardware kit

BASE WEIGHT ¹	79 lb (35.8 kg)
MAXIMUM SCOOTER/PC OVERALL WEIGHT	150 lb ² (68 kg)
MAXIMUM PRODUCT WIDTH	24" (61 cm)
OVERALL PLATFORM LENGTH	39" (99.1 cm) 37" (94 cm) max. product length
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21" (53.3 cm)
MOTOR	Fully sealed 12-volt DC
WARRANTY	3-Year Transferable
HITCH CLASS	II, III, IV
OPTIONS	• Battery pack w/charger • Swing-away adaptor
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover

1. Does not include any accessories or hardware.

2. Lifting weight does not include accessories, upgraded seating and/or hardware.

