

# VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER

NOW  
AVAILABLE

TRUE-INFINITE  
POSITIONING AND  
MEMORY TOGGLE  
HAND CONTROL!



Shown in Badlands Walnut Fabric



COUNTLESS  
POSITION  
COMBINATIONS

## ELEGANCE 2¹

A truly elegant, sophisticated recliner is more than just physically alluring: It's the precise combination of function, style and comfort. The **Elegance 2 Collection** from VivaLift!® Power Recliners is all of these and more, providing a cozy place to relax with your family.

The Elegance 2 has been upgraded with an ergonomically designed memory toggle remote that features two programmable memory position buttons along with up/down/recline functions for customizable comfort.



Scan here for  
more information  
on **VivaLift!®**









**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

# ELEGANCE 2<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION<sup>1</sup>

COLLECTION	ELEGANCE 2 <sup>1</sup>	ELEGANCE 2 <sup>1</sup>
MODEL	PLR-975M	PLR-975L
POSITION	True-Infinite	True-Infinite
HEIGHT RANGE	5'4"-5'9"	5'10" and above
SIZE	Medium	Large
WEIGHT CAPACITY	400 lbs.	400 lbs.
SEAT-TO-FLOOR	20"	20.5"
SEAT DEPTH	21"	22"
SEAT WIDTH	21"	22"
TOP OF BACK TO SEAT	27"	29"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	43"
DISTANCE FROM WALL	14"	16"
OVERALL WIDTH	35.5"	37.5"
CHAIR WEIGHT	152 lbs.	157 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated


## INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

## STANDARD FEATURES

- NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- Lithium battery backup ensures your recliner is operational when the power goes out.

## STANDARD FABRIC: BADLANDS



MUSHROOM

WALNUT

STEEL

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**ES QUICK SHIP** = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>