




AmeriGlide



platinum

STAIRLIFTS



Pictured
Inside Curved Rail

- When folded, the stairlift measures 17.3":
The diameter of the Platinum is unbeatable
- It is always ready to use, thanks to the battery
and the automatic Sleepmode system

The Platinum is the most compact.



platinum

STAIRLIFTS

Since 1996 the Platinum stairlift has been a leader of style and quality.

Available in:

- Canada
- United States
- Argentina
- United Kingdom
- France
- Italy
- Spain
- Portugal
- Austria
- Switzerland
- Germany
- Denmark
- Holland
- Poland
- Estonia
- Finland
- Sweden
- Greece
- Turkey
- Saudi Arabia
- South Africa
- India
- Australia


- THE CUSTOM MADE DOUBLE ROLLER AND TWIN RAIL SYSTEM GUARANTEES MAXIMUM STABILITY ALONG THE TRACK
- THE ARMS, SEAT, AND FOOTREST CAN BE FOLDED UP AND DOWN INDEPENDENTLY TO SUIT YOUR NEEDS.

Track supports

Attention is paid to every detail during manufacturing of a Platinum Stairlift.

Pictured

Outside Curved Rail

A man in a light blue shirt and dark trousers is sitting on a compact, white stairlift on a staircase. The stairlift is mounted on the wall and has a seat and footrest. The man is smiling and looking to the right. The staircase has wooden steps and a white handrail. The walls are wood-paneled. A potted plant is visible in the bottom right corner.

“The stairlift I chose for my house is convenient for everyday life. The ergonomically designed seat and its adjustment are truly innovative. This is a stairlift without equal.”

We chose the Platinum because:
it is the most compact stairlift in the world
it is the only one entirely adjustable for my needs
it is the only one with a diagnostic system built into the armrest
it is the only one with removable and washable seat covers.



“My father is an active person, the stairlift helps him to use the stairs and thus continue to enjoy a full life, allowing him to get on with his work and hobbies. I thought that all stair lifts would be the same, but no, I found that Platinum Stairlifts have benefits that other stairlift companies do not.”

for stairlifts choose Platinum.

This is a true story and location, Actors are used for the photograph.

Live life your way

Load

up to 352 lbs (264 lbs on standard model). The double roller and twin rail is the most compact in the world, it is secure and stable even with the maximum load.

Always safe

All Platinum stairlifts work with a battery, and the inVISION diagnostic display is easy to see so you know the status of your stairlift at all times.

Can be installed anywhere.

The Platinum stairlift is the ideal solution for any type of stairs, straight or curved, inside or outside the home.

In closed position

the Platinum stairlift takes up a minimum space. It can be parked at the foot or at the top of the stairs.



without worrying about the stairs.



Easy to use
thanks to the ergonomically designed controls the Platinum stairlift is user friendly.



Reaching the top
of the stairs the Platinum stairlift seat swivels to allow you to alight safely onto the landing.

Freedom

I re-discovered
the pleasure of freedom
and accessibility in my
home as the stairs were
becoming more difficult
all the time.

I am happy
because I have chosen a
stairlift with double rollers
and twin rail construction
which takes up the minimum
space on the stairs.
Platinum is the most
compact design in the world.

I am free
to leave the seat at the
bottom or top of the stairs
where I can call it using one
of the twin remote controls.

I move safely
because the battery is
always charged the
stairlift will still operate
even if there is a power
cut.

The security sensors
stop the stairlift safely if
there is an obstruction in
the way.



I prefer the Platinum.



The stairs, although they are not very wide they remain free and my children and grandchildren can use the stairs without hinderance.



With a curve at the bottom of the stairs the chair can be parked in the space available without blocking the doorway.

Tailor Made Track

The Platinum is designed for your stairs with attention to any obstructions to obtain optimal functionality and safety.



The Least Intrusive Track in the World

Before buying a stairlift check carefully how much space the track will take up, at the foot of the stairs. On the stairway, on the curve and on the landing. The Platinum can not be beaten.

The Key Features

Easy to use Joystick

A comfortable, natural and easy to use joystick. Even for those with limited mobility.

Seat covers removable and washable to keep your stairlift clean and tidy.



Exclusive to Platinum

Exclusive to Platinum

InVISION display set in the armrest so it's easy to see and understand the functions of your stairlift.

Exclusive to Platinum

Safety harness for peace of mind.

Swivel seat reaching the top of the stairs you can alight straight onto the landing safely.

Emergency hand winder In an emergency the seat can be moved using a manual handle.

ON/OFF switch to switch the stairlift on and off as required.

12 Security sensors any one sensor will stop the chair safely if an obstruction is detected.

Isolating Key located under the armrest.

Footrest lever to lift and lower the footrest without having to bend down.

Exclusive to Platinum

ERGNOMICALLY designed seat studied to optimise the position of the legs whilst seated.

unique in the world.



1

The only stairlift with light touch joystick control designed to fit the natural position of the fingers.



2

The only stairlift with the InVISION display in the armrest. In case of technical problems the display identifies the precise problem with a code number, enabling a swift resolution to the problem.




3

The only stairlift with ergonomical adjustment
The Platinum is adjustable for your comfort. It has 12 separate adjustments that be made to the seat and armrests. This is a stairlift without equal.



4

The only stairlift with completely removable seat covers
It's so easy to remove and wash the covers which is a necessity and not a luxury. This why we include this feature for the Platinum stairlift.



To lower the footrest is simple and it adjusts independently of the chair arms.

The two remote controls supplied with every stairlift allow you to call or move the seat remotely.



- Continues to work even in a power outage
- Can carry up to 264 lbs. (352 lbs - HD Model)

Exclusive to Platinum

- Start delay to move off safely after a few moments
- The speed of the curve is set independently.

Moving up and down the stairs

in total safety



The safety key under the **armrest** disables the stairlift to avoid accidental start by children.



The security sensors set in the seat, the footrest and the carriage will stop the stairlift safely if an obstacle is detected.



also outside the home.



CLASSIC seat

with all Platinum stairlifts.
Quality fabric with a tailored finish.
When you purchase a platinum you will always receive two sets of seat covers



ERGO seat. A patented design

a more compact seat with more natural posture for the legs. A lumbar support option is also available.

With a tailor made stairlift

according to my needs and preference.

Fabric - Red



You can choose the fabric and color that you prefer.

Fabric - Blue



Vinyl - Dark



Fabric - Orange



Vinyl - Sand




Fabric - Grey



Custom - (optional)

By request you can ask for the fabric which you like best. Truly custom made.

**Exclusive
to
Platinum**



Every Platinum Stairlift is individual



A precise study, attention to detail is essential to build a stairlift custom made for you.

DETAILED STUDY OF YOUR STAIRWAY



CUSTOM MADE

How we can help you choose the best option for your home?

Call Us

to request a free assessment by our consultant in your home.

A personal consultation

to explain the benefits of each product and find the best solution for you.

We install

the stairlift in less than a day, without any modification to the stairway.

We come to your home

quickly to evaluate your specific needs and solutions.

Custom made

your stairlift will be designed to perfectly fit your staircase.

Your stairlift

is ready for use immediately to make your day to day life easier.



Showroom



COMPUTER CONTROLLED MANUFACTURE



Attention to detail
The guide rollers and track are built and tested by specialists with great attention to detail.

ERECTED AND TESTED IN THE FACTORY BEFORE DELIVERY

Guarantee and aftersales service

Aftersale service

An efficient after sales service is a fundamental element in the life of a stairlift.

Quick response

Due to the network of AmeriGlide dealers throughout the country clients can count on quick professional service.

Diagnostic DISPLAY

Thanks to the InVISION diagnostic display on all AmeriGlide stairlifts our technicians can also resolve problems with a telephone call, reducing your waiting time.

Guarantee

covered by a 10-year warranty against defects in the drive train and its component parts, and a 5-year warranty against defects in all other stair lift parts and components.



Technical Assistance



For all these reasons. Today I have a Platinum stairlift in my home.

