

Deluxe Flame Rollator with 8-inch Wheels

The **ProBasics Deluxe Flame Rollator** represents both flair & function for those who need a little more personality from their walking aid. Large 8" wheels offer superior maneuverability and are ideal for both indoor and outdoor use. Constructed of lightweight aluminum and equipped with both a removable flip-back padded backrest and built-in strap for easy foldability, this rollator is designed to accompany and take your patient further whenever on the move. 300 lb weight capacity.







Water-Resistant Easy-to-Operate Nylon Pouch Ergonomic Hand Brakes

8-inch Wheels

FEATURES

- Lightweight aluminum frame in burgundy and blue flame finish
- Two-button removable flip-back curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable handles for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with soft grip tires, ideal for indoor and outdoor use
- HCPCS* Code: E0143

E0156 (seat)

COLORS

Burgundy Flame Finish



SPECIFICATIONS

| SILCIIICATIONS | | | |
|------------------------------|--|--|--|
| Overall Width: | 24" | | |
| Width b/w Hand Grips: | 18.25" | | |
| Seat Dimensions: | 14" (W) × 12" (D) | | |
| Seat-to-Floor Height: | 20.5" | | |
| Weight Capacity: | 300 lb | | |
| Adjustable Handle Height: | 31" – 36" | | |
| Wheel Size: | 8" | | |
| Product Weight: | 15.4 lb | | |
| Carton Dimensions: | 24.8" × 9.1" × 35" | | |
| Shipping Weight: | 18.3 lb | | |
| Warranty: | Limited lifetime on frame One year on all parts | | |
| | | | |

| ltem # | Description | Color | UPC | UOM |
|---------|--|----------------|--------------|------|
| RLA8BGF | ProBasics Deluxe Rollator with 8-inch Wheels | Burgundy Flame | 815067072521 | 1/cs |
| RLA8BLF | ProBasics Deluxe Rollator with 8-inch Wheels | Blue Flame | 815067072538 | 1/cs |

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.